



► CUBIC 18

◀ کویک ۱۸

CUBIC 18 Liters Class B Autoclave	
Dimensions	W540×H455×D660 mm
Chamber Size	W220×H220×D400 mm
Maximum working pressure	2.25bar
Minimum vacuum	-0.8bar
Maximum working temperature	138°C
Hydrostatic pressure	4bar
Dimension of trays	6pcs W205×H14×D365 mm
Working temperature	134°C~121°C
Vacuum pump	Diaphragm
Power supply	AC220 V 50-60Hz
Power	2100W
Control system	Microprocessor
Printer	Thermal

