



► CUBIC 25

◀ کویک ۲۵

CUBIC 25 Liters Class B Autoclave	
Dimensions	W560×H455×D660mm
Chamber Size	W250×H250×D400mm
Maximum working pressure	2.25bar
Minimum vacuum	-0.8bar
Maximum working temperature	13۵°C
Hydrostatic pressure	4bar
Dimension of trays	7pcs, W240×H14×D365 mm
Working temperature	134°C~121°C
Vacuum pump	Diaphragm
Power supply	AC220 V 50-60Hz
Power	2100W
Control system	Microprocessor
Printer	Thermal

