



MOTOR (ISM -B70L)

Model	ISM-B70L
Max.Speed	40,000 RPM
Max.Torque	70 N.cm(20:1)
Internal Resistance	730mΩ
Internal Inductance	100uH
Dimension	24mm(W) x 100mm(W)
Weight	156g
Optic	Yes

CONTROLLER (ISE-270M)

Model	ISE-270M
Power Supply voltage	AC110 ~ 240V
Frequency	50/60Hz
Power Consumption	150VA
Max. Coolant flow rate	110ml/min
Fuse	2 x 250V - T2.0AH
Dimension	285mm(W) x 250mm(W) x 120mm(H)

FOOT PEDAL (FC-V2)

Model	FC-V2
Speed control	Variable
Foot control Function	
	Program
	Coolant
	Clockwise/Counter clockwise