

DTE[®]

The choice of millions of dentists worldwide

GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.

Registration Add & Manufactory Add:
Information Industrial Park, Guilin National High-Tech
Information Park 541004 Guilin, Guangxi China
Tel:
Europe Sales Dept.: +86-773-5873196, +86-773-2125222
North America, South America &
Oceania Sales Dept.: +86-773-5873198, +86-773-2125123
Asia & Africa Sales Dept.: +86-773-5855350, +86-773-2125896
Fax: +86-773-5822450
E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com
Website: <http://www.glwoodpecker.com>

Distributor:

Scan and Login website
for more information



DTE[®]

The choice of millions
of dentists worldwide

Implant-X

Precision strength safety

Implant-X Your favorite assistant of high-precision
dental implanting treatment



Outstanding quality makes it your solid choice

Implant-X adopts mute motor imported from Swiss.
Advanced FOC controlling technique blends low speed performance with high speed performance.
In order to verify the clinical features of flexible response and stable safety,
Implant-X R & D personnel go into the operating room,
getting them involved in clinical applications.
Under the instruction and assistance of experts from home and abroad,
they finally complete the clinical optimization of product.
DTE spares no efforts to provide you high quality products.



Intuitive user interface

Wide color screen with simple and delicate interface
Eight modes for option; accommodate all implanting system
Color touch screen and glass surface bring you better operating experience
Accurate torque control and real-time display ensure surgical safety

Motor with precision and strength

Imported powerful Swiss motor

6.2Ncm motor torque can guarantee 80Ncm output torque,
making difficult surgery easier

Motor and cable can be autoclaved.

Lighter micro motor realize flexible operation

LED motor enables more efficient and accurate treatment.



Strong and stable contra-angle

- LED contra-angle ensures visibility, helping to achieve treatment of high accuracy and efficiency.
- High precision contra-angle with tiny oscillation enables stable torque output under high speed and low speed.
- Exported advanced bearings guarantee performance and durability of contra-angle.



Color touch screen with glass surface

The colour touch screen with its glass surface makes operation of the control unit simple and intuitive. The glass surface is easy to clean. Wiping disinfection ensures safety in treatment area.

Flexible multi-function foot pedal

Multi-function foot pedal can realize aseptic operation and liberate two hands.



Mute peristaltic pump

Mute peristaltic pump water supply system is easy to connect. Be compatible with disposable oral water pipe to ensure safety of treatment.



Perfect balance

Balancing between weight of contra-angle and motor guarantees better holding

Perform prolonged operation without fatigue

LED contra-angle ensures visibility, helping to achieve treatment of high accuracy and efficiency.





DTE

SPEED 26000

TORQUE 50

PROGRAM 1

SPRINT
TORQUE
PROGRAM
LIGHT
MEMORY
COOLANT
EMERGENCY