

Detector	CMOS
Pixel Size	29.6 μm
Theoretical Resolution	17 lp/mm
Dynamic Range	12 bit
Active Area(WxL)	Size 1.0 _ 20 x 30 mm Size 1.5 _ 24 x 33 mm Size 2.0 _ 26 x 36 mm
Dimensions(WxLxT)	Size 1.0 _ 25.4 x 36.8mm (1.00 x 1.45 inch) Size 1.5 _ 29.2 x 39.5mm (1.14 x 1.55 inch) Size 2.0 _ 31.3 x 42.9mm (1.23 x 1.69 inch)
Thickness	4.8 mm(0.19 inch)
Cable Length	2.7 m

