

Classic LED curing light

NEW BLUEDENT POWER PEN CORDLESS

| | |
|------------------------------|--|
| Wavelength range | 450 - 490 nm |
| Curing mode - Timer | RAMP - 10/20/30 sec. Sound signal every 10 sec. |
| Light intensity | up to 1000 mW/sq.cm |
| Weight / Length of handpiece | 100 g / 145 mm |
| Battery | Li-Ion 3.7 V, 800 mAh |



BLUEDENT LED POWER PEN



| | |
|------------------------------|---|
| Wavelength range | 450 - 490 nm |
| Curing mode - Timer | RAMP - 10/20/30 sec. Sound signal every 10 sec |
| Light intensity | up to 1300 mW/sq.cm |
| Weight / Length of handpiece | 80 g / 142 mm |
| Power supply | 220VAC / 24VAC / 12VAC |

 **BG LIGHT LTD** BULGARIA
www.bglight.com