

Polymerisation + bleaching in 1 product

2 in 1


BLUEDENT LED SMART CORDLESS

Wavelength range	450 - 490 nm
Modes	RAMP, FAST, PULSE, LONG
Timer	10/20/30 sec. / 9 min.
Light intensity	up to 1300 mW/sq. cm
Weight / Length of handpiece	230 g / 185 mm
Weight of station	540 g
Battery	Li-Ion 3.7 V, 2400 mAh (up to 4000 sec. work)
Charging base with 4-level light meter	Option



BLUEDENT LED SMART cable (built-in) – LED curing light

Wavelength range	450 - 490 nm
Modes	RAMP, FAST, PULSE, LONG
Timer	10/20/30 sec. / 9 min.
Light intensity	up to 1300 mW/sq. cm
Weight / Length of handpiece	230 g / 185 mm
Power supply	220V AC / 24V AC / 12V AC



BLUEDENT LED PRO

Wavelength range	450 - 490 nm
Curing mode - Timer	RAMP - for polymerisation: 10/20/30 sec. LONG - for continuous work during teeth bleaching and brackets fixing: 30 min.
Light intensity	up to 1400 mW/sq. cm
Weight / Length of handpiece	120 g / 185 mm
Power supply	220V AC / 24V AC / 12V AC

