
# Multiple shooting methods







Remote control with a tripod and extended exposure hand brake line

# Technical Parameter Portable X-ray machine

			_
Tube Voltage	Highest Voltage	Tube Current	
60KV	70KV	2mA	
Exposure Time	Focus	Anode Target Angle	
0.01-2s	0.4mm	12°	
Total Filtration	Battery Capacity		
1.0mmAl	950mAh		

#### Fussen Technology Co.,Ltd


- 2 1F,3#Building ,Tongfuyu Industrial Zone, Taoyuan Street,Nanshan District, Shenzhen,China
- **(**) Tel:+86 755 2167 0890 813
- $\ \, \bigoplus \, \mathsf{Website} \colon \mathsf{www.fussengroup.com}$

# **FUSSEN**



Portable X-ray

A smart dental digital solution



#### High frequency DC technology

Ultra-low radiation, the exposure dose is reduced by more than 30% compared with the standard generator, maximizing effective rays and greatly reducing the radiation intensity




# High capacity lithium battery

950 mAh ultra high capacity battery, battery can be replaced, low cost of use



# Humanized and easy operation

Shoot different teeth, fine-tune the shooting time, and achieve personalized and precise control



# Small focus to get high definition images

High-voltage generator with 0.4mm focus, high-definition imaging, enjoy top quality

