



# UltraMint & UltraMint Pro

All-round Ultrasonic Scaler

## All-round Ultrasonic Scaler

Scaling | Periodontic | Endodontic



### Soft Perio Power for Ultimate Comfort

Autoclavable plastic box - tips\*6, universal wrench\*1, EMS type handpiece\*1

Autoclavable plastic box - tips\*6, universal wrench\*1, Satelec type handpiece\*1



\*Applicable for all tips from Eighteenth


Smart touch system-ultrasensitive & waterproof



Intelligent self-cleaning-makes disinfection easy and avoids cross-infection





Two ways of water supply- Changing supply methods is easier than ever before



900ml stand-alone bottle with clear 20ml scale mark makes it convenient for formulating liquid.




Detachable handpiece with LED- autoclavable, ergonomic, brighter visibility



Detachable and autoclavable handpiece holder, resists cross-infection



Water control directly from handpiece



Ergonomic foot pedal



• EMS type

• Satelec type

## One scaler, 2 solutions

Wide range of frequency: 25kHz-42kHz

– We are proud to give you two solutions to perform at your best—both handpieces (EMS type and Satelec type) provided by Eighteeth are compatible with UltraMint & UltraMint Pro.\*  
Scaling can be as convenient as you can imagine.

\*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.

## AI Engine — Scaling has never been so intelligent!



Customized treatment based on tooth surface conditions



Smart touch control system—ultrasensitive & waterproof



Elliptical stroke patterns with polishing effect




Intelligent degenerative feedback



Digital controlled automatic frequency tracking  
Ensures machine works on optimum frequency



Sine wave drive  
Avoids simulation caused by instantaneous current output; ensures comfortable treatment




Tips: G1 G2 G4 P1 P3 E4 universal wrench, EMS type handpiece



Tips: E3-SD E5-S E7-S P4-S G1-S G2-S universal wrench, Satelec type handpiece

## Scaling



## Periodontic



## Endodontic



## Endodontic



## Cavity Preparation



\*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.

\*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.