



HyperLight

Portable X-ray unit




HyperLight

Engineered for accurate images

- Clinical Imagery



- Sharp, high resolution images with a **0.4mm** focal spot



Designed with safety in mind



Built-in leaded double shield to protect the operator from leakage radiation



Locking system to prevent unintentional exposure



Specifications

Model	HyperLight	Weight	1.8 kg
Focal Spot	0.4 mm	Battery	2500 mAh x 4
Tube Voltage	65 kV	Exposure Time Range	0.02-2 Sec
Tube Current	2.5 mA	Line voltage	100V-240V, 50Hz/60Hz

HyperLight




HyperLight



Simplify your digital workflow

- Compact size and lightweight can make you operate it by one hand only
- Intuitive user interface
- Smart angulation system



Maximize your diagnostic potential

Perfect combination



HyperLight