



MEDIT *i900* classic

Effortless Precision, Your Way

MEDIT i900 classic

Introducing the Medit i900 classic—a refined, advanced scanning solution built from your insights. With a reimagined tactile button interface that offers intuitive, quick, and reliable operation, it transforms the scanning process into an effortless experience. Designed to enhance clinical workflows, the Medit i900 classic combines cutting-edge technology with user-focused innovation, ensuring precision and flexibility tailored to your unique needs.

MEDIT i900



MEDIT i900 classic



Effortless Precision,

Your Way



MEDIT i900 classic



MEDIT i900

Advanced or classic?

MEDIT

MEDIT *i900 classic*

**PRECISION IN
EVERY CLICK**

Simpler operation

- Intuitive button interface
- Quick, reliable response



The only touch scanner

- 360° Touch Band & Touch Pad
- Just lift and start scanning

MEDIT *i900*

**SEAMLESS
TOUCH EXPERIENCE**

MEDIT Intraoral Scanner Lineup

MEDIT

Introducing the MEDIT i900 classic to our scanner lineup

MEDIT *i*900 | MEDIT *i*900 classic


FLAGSHIP

MEDIT *i*700 | MEDIT *i*700 wireless

PERFORMANCE

MEDIT *i*600

ESSENTIAL



The MEDIT i900 Family

Effortless Precision, Your Way



The Medit i900 Family

Two ways to scan. One seamless experience.

MEDIT i900

Touch



MEDIT i900 classic

Button

One seamless unibody, hygienic by design

Experience unmatched durability and elegance with the Medit i900 classic's unibody construction. By eliminating assembly seams, it remains pristine and resists dirt and disinfectants. Engineered to reduce cross-contamination, its seamless design not only prioritizes safety but also delivers perfect weight distribution for effortless handling.



Compact & lightweight

MEDIT

Plug-and-play scanning, scan longer with ease.

Designed for effortless mobility, the MEDIT i900 classic weighs just 165g—lighter than your smartphone—and feels almost weightless in your hand. With its plug-and-play design and newly enhanced, ultra-reliable cable, this scanner delivers seamless performance every time. Embrace a scanning experience where power meets unparalleled freedom, transforming every session into a breeze.



New 3rd-gen optical engine

MEDIT

The new 3rd-gen optical engine drives our dual optical camera system, delivering unprecedented clarity and precision. This breakthrough technology enhances metal scanning performance and captures richer textures, ensuring every scan is detailed and accurate. With deeper scanning depth and improved image quality, our advanced optical solution redefines digital dentistry, setting a new benchmark in precision and performance.



MEDIT *i900* | MEDIT *i900 classic*


MEDIT *i700*

Dual optical cameras

MEDIT

Double the clarity, twice the precision.

Challenge reflective metals and navigate the varied environments of the oral cavity with ease. Our advanced technology ensures enhanced image quality and accuracy, turning once-difficult scans into seamless, precise results every time.




With the capability to scan in tighter and deeper spaces with remarkable accuracy, our scanner provides unparalleled performance, ensuring that every detail is captured with the utmost clarity.

Capture more, faster


MEDIT

Embracing a larger scanning area, the MEDIT i900 classic efficiently captures extensive data, reducing both scan time and processing duration. Unlike other products that may suffer data interruptions based on the scan path, the MEDIT i900 operates smoothly and rapidly—even with challenging edentulous models. It transcends conventional limitations, ensuring every scan is continuous, coherent, and a testament to exceptional performance.



MEDIT i900 | MEDIT i900 classic

18 x 15 mm (552 x 460 px)



MEDIT i700


Wider field of view

MEDIT

Capture more in a single sweep.


Experience effortless, comprehensive scanning with our enhanced field of view.

Capture a wider area with ease, ensuring speed and precision in every session.



MEDIT i900 | MEDIT i900 classic

18 x 15 mm (552 x 460 px)




MEDIT i700

Easy for beginners with AI Scan Assist

MEDIT

Effortless Scanning with the Reliability Map Feature

Our innovative Reliability Map feature provides real-time precision checks during scanning. It identifies missing spots and displays arrows to guide you. Optimal scan areas turn greener, ensuring reliable results every time.



Not just accessories

3 Types of Tips Available

Medium Tips

A versatile choice, offering the perfect balance of coverage and control for everyday scanning.



Small Tips

Designed for precision in tight spaces, ensuring comfort and accuracy when scanning smaller arches or hard-to-reach areas.



Large Tips

Capture more in less time—maximize scan efficiency and scan edentulous spaces with ease.



Dual optical cameras

- New 3rd-gen optical engine
- Improved metal scanning
- Enhanced texture quality



One seamless unibody

- Hygienic and easy to clean
- Durable and innovative interface
- Enhanced response performance



Tactile Button

Precision in Every Click



Touch Band & Touch Pad

Seamless Touch Experience

Compact & lightweight

- Plug-and-play scanning
- Lighter than a smartphone
- Ultra-reliable cable




Wider field of view

- Faster scans
- Fewer image captures
- 18 x 15mm FOV size



Smart scan assistant

- 30mm deep scanning
- AI-powered scan guide
- Fast preheating system



MEDIT Shield Program



KEEP CALM AND SHIELD ON

Protect your investment with Medit Shield Program.

About Medit Shield Program

► **Unlimited Replacements:**

Accidental drop? No worries—
if your scanner takes a fall, we replace it.

► **Rapid, Direct Service:**

Count on fast, reliable support whenever you need it.

► **Comprehensive Coverage:**

Enjoy extensive protection and peace of mind
for your investment.



Medit SmartX

Effortless All-on-X Workflow

- **No extra scanner needed**—
fully compatible with i900 Family, i700, and i600
- Smarter and faster All-on-X workflow
- Optimized alignment for precision
- Real-time library matching
- Scan body system continuously adapts
- Customizable workflow tailored to user




Open up a world of possibilities,



**ORTHODONTIC
SOLUTION**



**CHAIRSIDE
SOLUTION**



**IMPLANT
SOLUTION**



**DIAGNOSIS &
CONSULTATION**

across various dental solutions.

Infinite Possibility with MEDIT

At MEDIT, we're not just about creating products; we're about reimagining possibilities.

Our commitment to innovation and customer care propels the digital dentistry revolution, offering advanced 3D digital solutions that transform dental clinics and labs. With our intraoral scanners, desktop scanners, and sophisticated platform and collaboration software, we're setting new standards in dental technology.





Thank You

MEDIT