

Filter-FREE Innovation

Discover the **NEW DWX-43W** Wet Milling Solution



Roland DGA's DGSHAPE is the leading choice for lab-quality, esthetic crowns, bridges, onlays and veneers NOW introduces the **NEW DWX-43W** wet milling solution. With the latest filter-free innovation, combined with faster milling times for traditional ceramics and custom titanium abutments, your dental business will now produce superior same-day restorations and reduce operating costs while increasing ROI.

- **Open Architecture Eco-System Platform**
- **Filter-Free Technology**
- **Easy Integration with Ethernet LAN or USB Connectivity**
- **Designed Specifically for Wet Milling of Harder Materials Like: Chairside Zirconia & Amber Mill Direct, Powered by Roland DGA**

NEW DWX-43W Comes With PERKS

- Say YES to NO filters... **NEW DWX-43W** Filter-FREE Technology eliminates the need for coolant tank filters reducing the cost of use by 20%.
- Cuts minimal routine maintenance by 50% with the new automatic coolant jet optimization and draining cycle.

Added Value Extras:

- **NEW DGSHAPE** CAM software included with purchase.
- DGSHAPE authentic tools and accessories included with DWX-43W for expanded milling capabilities and greater ROI.

Faster Milling times, Increased efficiency.
The *NEW DWX-43W* is All-Milling-All-Gain!

For pricing and full details, visit www.rolanddga.com/DWX-43W

COMPATIBLE MATERIALS	Traditional & Pre-Crystalized Lithium Disilicates, Glass Ceramics, PMMA, Hybrid Ceramics, Chairside Zirconia-Powered by Roland DGA, Composite Resins, and Pre-Mill Titanium Blanks	
SPINDLE SPEED	15,000~60,000 RPM	
ROTARY AXIS	A	±360°
AIR PRESSURE REQUIREMENTS	0.2 Mpa or More	
AIR CAPACITY	39 L/min or more	
EXTERNAL DIMENSIONS	Width, Depth, Height	19 in. x 22 in. x 23.6 in.
WEIGHT	134.48 lbs (61 kg)	
ABUTMENT MATERIAL ADAPTER COMPATIBILITY	AK-1-ME (MEDENTIKA®), AK-1-DE (DESS®), and AK-1-GE (GeoMEDI®)	
AUTOMATIC TOOL CHANGER	Houses 6 Milling Tools	
ADDITIONAL FEATURES	Filter-less Coolant tank technology enclosure system and Automatic Draining Cycle	
ACCESSORIES INCLUDED	DGSHAPE CAM Software, 6-Pin Material Adapter (#ZV-42W), ZGB2-125D (3), ZGB2-125DS (1), ZGB2-75D (1), ZGB2-50D (3), ZGB2-25D (1), ZGB-125D (3), and ZGB-50D (3).	
INTERFACE	USB / LAN Connectivity	
LOADABLE WORKPIECE	3 Standard Mandrel Material CAD/CAM Blocks with Internal Adapter	
	6 Standard Mandrel Material CAD/CAM Blocks with included 6-Pin Adapter	

Warranty and Support

A two-year warranty plus support from a network of dedicated dealers and certified service engineers ensures that your machine is always up-and-running. From our dedicated customer support team to a wealth of training and education, DGSHAPE by Roland cares about your customer needs before and after you purchase your device.

Additional Features

Open-Architecture Platform: Allows users the flexibility to utilize their choice of restorative materials, scanners, and CAD Design Software. Users can customize a workflow that is truly built to build their business.

DGSHAPE's proprietary VPanel and DGSHAPE Cloud allow for monitoring from anywhere, anytime on any device.

Discover growth in your dental business with the new advanced milling strategies today. For more information contact your dedicated Roland DGA sales representative.

