

# AutoScan-DS-EX Pro(H)

Embrace The 5.0 MP Era



SHINING  
DEN



# Technical Specifications

Model	AutoScan-DS-EX Pro(H)
Camera Resolution	2*5.0 MP
Accuracy	≤ 8 μm
Scan Time	Bite scan: 7s; Upper/Lower jaw scan: 14s; 1~4 dies scan: 18s; 5~8 dies scan: 30s; Impression scan: 66s
Articulator Scan	Support articulator transfer
Output Format	STL, PLY, OBJ
Light Source	Blue Light
Dimension	385*260*270 mm
Weight	5.0 Kg
Interface	USB 3.0
Power Supply	DC 24V

Notice: SHINING 3D reserves the right to introduce modifications or alterations to any of the specifications and images used in this document.

Scanned with CamScanner