


Intraoral Disease Detector	Intraoral Camera		Curing Light	Apex Locator
-IC-WHCD100 -IC-DCDU75	-IC-WHT60 -IC-DCAM70 -IC-EZAM80 -IC-WHMA90	-IC-WHW61 -IC-DCAU71 -IC-EZAU81	-CL-DP40 -CL-DL10 -CL-DC20	-AL-DFA20 -AL-DF10



Macro to Full arch

Camera is able to have Macro to Full arch Image.



*EIOS When equipped,
it is possible to capture even
the face. (Extra Illumination Optic System)



3 Color Impressions

Possible to Select Color Impression(Original, Warm, Cool)



Software compatibility

Able to be used with almost any dental imaging software in the market.



GoodDoctors

Enjoy Wireless
Intraoral Camera



Whicam Air



Dr's cam Air

Air Direct

New conceptual Wireless Solution
for Medical and Dental environment



www.gooddrs.com
E-mail. info@gooddrs.com



Good Doctors Co., Ltd.
#208, B-dong, 283 Bupyeong-daero, Bupyeong-gu (Woolim Lions Valley,
Cheongchoen-dong), Incheon, 21315, Republic of Korea
Tel: +82 32 424 6325 Fax: +82 32 424 6326 GDBR-IC-AD-Air Direct_Rev.0

New conceptual Wireless Solution for Medical and Dental environment


A Simple selection and Easy connection

I Infinite Channel

R Real Plug and Play



Intraoral Camera



Chair 01 Chair 02 Chair 03 Chair 04 ... Chair

	Air Direct	UWB	Wifi
Connection Initiative	Camera to Receiver	Receiver to Camera	Receiver to Camera
Computer Operation when Connecting	Do not Need	Need	Need
No. of Channels	Unlimited	-	-
Driver Software	Not Required	Need	May Need
Management S/W Support	O	O	X
Transmission Distance	Long	Short	Long
Wireless Signal Penetration Power	Strong	Weak	Strong
PC Support	O	O	O
Mac Support	O	X	-



Whicam Air



Dr's cam Air

1 Features of Air Direct Intraoral Camera

- High speed digital wireless**
High speed digital wireless provide more clear live image without interruption.
- Macro to Full arch**
Camera is able to have Macro to Full arch Image.
- Auto focus**
Real Autofocus function shows clear images.
- HD image**
High definition show more detail images.
- A Simple selection and Easy connection**
Select the Connection signal on the camera.
- Infinite Channel**
Unlimited Channel Support.
- Real Plug and Play**
Not required Driver Software.
- Mac Support**
Mac and PC support with wireless.
- Wireless mouse**
Built in Gyro mouse.
- OLED screen**
Bright and Clear OLED screen indicates various device status.
- 280° twist head**
Easy to get approach anywhere in oral cavity.
- Adjustable LED brightness**
Can adjust LED brightness to avoid light reflecting.
- 3 Color Impressions**
Possible to Select Color Impression(Original, Warm, Cool)
- Software compatibility**
Able to be used with almost any dental imaging software in the market all over the world.
- Twain driver**
Twain Driver Provided.
- Big capacity battery**
Long time use with big capacity battery.
- Adjustable head angle (Whicam Air)**
User can use the camera with prefer angle.

2 Specification

Model name	IC-WHT60 (Whicam Air)	IC-DCAM70 (Dr's cam Air)
Camera	Operation voltage	3.7V
	Measurement	Approx. 266 x 31 x 37 mm
	Weight	Approx. 255 x 30 x 28 mm
	Resolution	Approx. 140g
	Sensor	5M pixels
Charger	Operation voltage	5V
	Measurement	Approx. 169 x 148 x 89 mm
	Weight	Approx. 110 x 105 x 68 mm
Disposable cover	Dimensions	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

3 Minimum System Requirement

PC system	3.3GHz CPU, 4GB RAM
OS	Windows 7
USB port	USB 2.0

* Upper specification is changeable without prior notice.