



Runyes Future Product Plan

V3 Pro / V5 / V6

For internal use only



V3 Pro



V5



V6

For internal use only

Specifications		V3 Pro	V5	V6
Handpiece	Dimession(L x W x H)	240x49.5x30.5mm	240x40x32.5mm	265x40x32.5mm
	Weight	210g	160g	220g
Tip	Standard Tip Size	79x19.7x15.8mm	115x40x32.5mm	115x40x32.5mm
	Mini Tip Size	80x16x11mm	115x40x32.5mm	115x40x32.5mm
	FOV(Standard)	15x15mm	16x14mm	16x14mm
	FOV(Mini)	12x12mm	12x12mm	12x12mm
	Depth of Field	15mm	25mm	25mm
	Autoclavable	100 times autocalvable		

For internal use only

Specifications		V3 Pro	V5	V6
Accuracy	Crown	10μm	10μm	10μm
	Full Arch	20μm	20μm	20μm
Scan Speed	Single Arch	30s	30s	30s
	Full Arch with Occlusion	90s	90s	90s
Calibration	Calibration Free			
Output	STL, OBJ, PLY			
Connectivity	Type C (USB 3.0 with connector)			
Languages	Chinese, English, French, Spanish, Russian, Portuguese, Italian, Greek, Polish, Romanian, Japanese, Korean, Turkish, Czech, German, Hungarian			

For internal use only

Computer Configuration

Basic & Recommended Computer Specification Microsoft Windows 10/11, 64bit or above			
CPU	Desktop: Intel i7-9700 quad-core or above; (Recommended: Intel Core i7-11700) Laptop: Intel i7-9700H quad-core or above; (Recommended: Intel Core i7-11700H)	System	Windows 10/11 Pro/Corporate Edition
Graphics Card	NVIDIA GeForce 1660GTX or above, more than 6G (Recommended: NVIDIA GeForce 3060 or above graphics card with 6G or more memory) Note: AMD graphics cards are not supported.	Monitor Resolution	1920*1080
Memory	More than 16 GB	Port	USB 3.0 or above
Hard Disk	More than 256GB SSD or 128GB+1TB mechanical hard disk	Power Consumption	25 V/A

For internal use only



	V3 Pro	V5	V6
Connecting Method	One Cable	One Cable	Cable or Wireless

Wireless innovation for seamless dental scans


For internal use only




	V3 Pro	V5	V6
Depth of Field	15mm	25mm	25mm

Capturing every detail for superior care
Larger scanning area for a faster process

For internal use only



Heating >>> Fan




	V3 Pro	V5	V6
Anti-Fogging Method	Heating	Fan	Fan

For internal use only

Better cross infection avoidance

Longer & Cheaper ! !



	V3 Pro	V5	V6
Tip	Sterilization for 100 times	Sterilization for 100 times,Longer, Less cost	Sterilization for 100 times,Longer, Less cost

For internal use only



Ergonomic design

- Slimmer shape
- Better grip



	V3 Pro	V5	V6
Handle (Estimated Weight)	210g	160g	220g



Charge Using

- Battery charger
- Charging Dock(Type C Connection)



	V3 Pro	V5	V6
Recharging Method	X	X	Two ways of recharging



Battery Charging





Cable Connection Mode

